

Cat. #	品名	分子量	供給量	参考価格
45920	H-4-Abz-2-Chlorotrityl Resin	0.3~0.8 mmol/g	25g	30,400
	100~200 mesh, 1% DVB		100g	91,000
45913	Fmoc-Acp-2-Chlorotrityl Resin	0.3~0.8 mmol/g	25g	21,000
	100~200 mesh, 1% DVB		100g	58,300
44001	H-Ala-2-Chlorotrityl Resin	0.3~0.8 mmol/g	5g	8,200
	100~200 mesh, 1% DVB		25g	28,600
44007	H-Ala-2-Chlorotrityl Resin	0.3~0.8 mmol/g	5g	8,200
	200~400 mesh, 1% DVB		25g	28,600
44002	H-(N-Me)Ala-2-Chlorotrityl Resin	0.3~0.8 mmol/g	5g	35,000
	100~200 mesh, 1% DVB		25g	122,500
44003	H-β-Ala-2-Chlorotrityl Resin	0.3~0.8 mmol/g	5g	16,400
	100~200 mesh, 1% DVB		25g	57,200
45904	H-D-Ala-2-Chlorotrityl Resin	0.3~0.8 mmol/g	5g	8,200
	100~200 mesh, 1% DVB		25g	32,700
44004	Fmoc-Ala-2-Chlorotrityl Resin	0.3~0.8 mmol/g	5g	8,900
	100~200 mesh, 1% DVB		25g	30,900
44006	Fmoc-Ala-2-Chlorotrityl Resin	0.3~0.8 mmol/g	5g	8,900
	200~400 mesh, 1% DVB		25g	30,900
44005	Fmoc-β-Ala-2-Chlorotrityl Resin	0.3~0.8 mmol/g	25g	15,200
	100~200 mesh, 1% DVB		100g	52,500
44101	H-Arg(Pbf)-2-Chlorotrityl Resin	0.3~0.8 mmol/g	5g	14,400
	100~200 mesh, 1% DVB		25g	57,600
44102	H-D-Arg(Pbf)-2-Chlorotrityl Resin	0.3~0.8 mmol/g	5g	29,200
	100~200 mesh, 1% DVB		25g	105,000
44103	Fmoc-Arg(Pbf)-2-Chlorotrityl Resin	0.3~0.8 mmol/g	5g	29,200
	100~200 mesh, 1% DVB		25g	105,000
44201	H-Asn-2-Chlorotrityl Resin	0.3~0.8 mmol/g	5g	8,400
	100~200 mesh, 1% DVB		25g	33,600
44205	H-Asn-2-Chlorotrityl Resin	0.3~0.8 mmol/g	5g	8,400
	200~400 mesh, 1% DVB		25g	33,600
44202	H-Asn(Trt)-2-Chlorotrityl Resin	0.3~0.8 mmol/g	5g	9,600
	100~200 mesh, 1% DVB		25g	38,500
44203	Fmoc-Asn(Trt)-2-Chlorotrityl Resin	0.3~0.8 mmol/g	5g	29,200
	100~200 mesh, 1% DVB		25g	105,000
44204	Fmoc-Asn(Trt)-2-Chlorotrityl Resin	0.3~0.8 mmol/g	5g	29,200
	200~400 mesh, 1% DVB		25g	105,000
44302	H-Asp-OAll-2-chlorotrityl Resin	0.3~0.8 mmol/g	5g	7,000
	100~200 mesh, 1% DVB		25g	24,700
44301	H-Asp(OtBu)-2-Chlorotrityl Resin	0.3~0.8 mmol/g	5g	9,600
	100~200 mesh, 1% DVB		25g	38,500
48109	Fmoc-Asp(OAll)-2-Chlorotrityl Resin	0.3~0.8 mmol/g	5g	41,200
	100~200 mesh, 1% DVB		25g	72,000
44303	Fmoc-Asp(OtBu)-2-Chlorotrityl Resin	0.3~0.8 mmol/g	5g	11,900
	100~200 mesh, 1% DVB		25g	47,900

Cat. #	品名	分子量	供給量	参考価格
44304	Fmoc-Asp(OtBu)-2-Chlorotrityl Resin	0.3~0.8 mmol/g	5g	11,900
	200~400 mesh, 1% DVB		25g	47,900
45914	Fmoc-Cha-2-Chlorotrityl Resin	0.3~0.8 mmol/g	25g	17,500
	100~200 mesh, 1% DVB		100g	56,000
45918	H-Chg-2-Chlorotrityl Resin	0.3~0.8 mmol/g	5g	21,200
	100~200 mesh, 1% DVB		25g	31,600
44401	H-Cys(Acm)-2-Chlorotrityl Resin	0.3~0.8 mmol/g	5g	12,800
	100~200 mesh, 1% DVB		25g	51,100
44402	H-Cys(Trt)-2-Chlorotrityl Resin	0.3~0.8 mmol/g	5g	10,500
	100~200 mesh, 1% DVB		25g	36,800
44405	Boc-Cys(Acm)-2-Chlorotrityl Resin	0.3~0.8 mmol/g	5g	10,700
	100~200 mesh, 1% DVB		25g	26,300
44407	Fmoc-Cys(Trt)-2-Chlorotrityl Resin	0.3~0.8 mmol/g	5g	12,300
	200~400 mesh, 1% DVB		25g	22,900
44404	H-D-Cys(Trt)-2-Chlorotrityl Resin	0.3~0.8 mmol/g	5g	14,000
	100~200 mesh, 1% DVB		25g	56,000
44406	H-D-Cys(Trt)-2-Chlorotrityl Resin	0.3~0.8 mmol/g	5g	14,000
	200~400 mesh, 1% DVB		25g	56,000
44501	H-Gln-2-Chlorotrityl Resin	0.3~0.8 mmol/g	5g	8,400
	100~200 mesh, 1% DVB		25g	33,600
44502	H-Gln(Trt)-2-Chlorotrityl Resin	0.3~0.8 mmol/g	5g	9,600
	100~200 mesh, 1% DVB		25g	38,500
44503	Fmoc-Gln(Trt)-2-Chlorotrityl Resin	0.3~0.8 mmol/g	25g	17,500
	100~200 mesh, 1% DVB		100g	52,500
44601	H-Glu-2-Chlorotrityl Resin	0.3~0.8 mmol/g	5g	9,700
	100~200 mesh, 1% DVB		25g	31,800
44602	H-Glu(OtBu)-2-Chlorotrityl Resin	0.3~0.8 mmol/g	5g	9,600
	100~200 mesh, 1% DVB		25g	38,500
48108	Fmoc-Glu(OtBu)-2-Chlorotrityl Resin	0.3~0.8 mmol/g	5g	9,700
	100~200 mesh, 1% DVB		25g	29,200
48110	Fmoc-Glu(OtBu)-2-Chlorotrityl Resin	0.3~0.8 mmol/g	5g	9,700
	200~400 mesh, 1% DVB		25g	29,200
44701	H-Gly-2-Chlorotrityl Resin	0.3~0.8 mmol/g	5g	9,700
	100~200 mesh, 1% DVB		25g	32,100
44705	H-Gly-2-Chlorotrityl Resin	0.3~0.8 mmol/g	5g	9,700
	200~400 mesh, 1% DVB		25g	32,100
44703	H-Glycinol-2-Chlorotrityl Resin	0.3~0.8 mmol/g	5g	29,200
	100~200 mesh, 1% DVB		25g	105,000
44706	H-Glycinol-2-Chlorotrityl Resin	0.3~0.8 mmol/g	5g	29,200
	200~400 mesh, 1% DVB		25g	105,000
44704	Fmoc-Gly-2-Chlorotrityl Resin	0.3~0.8 mmol/g	5g	9,700
	100~200 mesh, 1% DVB		25g	23,400
44707	Fmoc-Gly-2-Chlorotrityl Resin	0.3~0.8 mmol/g	5g	9,700
	200~400 mesh, 1% DVB		25g	23,400

Cat. #	品名	分子量	供給量	参考価格
44801	H-His(Trt)-2-Chlorotrityl Resin	0.3~0.8 mmol/g	5g	9,600
	100~200 mesh, 1% DVB		25g	33,500
44803	Fmoc-His(Trt)-2-Chlorotrityl Resin	0.3~0.8 mmol/g	5g	11,700
	100~200 mesh, 1% DVB		25g	40,900
44802	H-D-His(Trt)-2-Chlorotrityl Resin	0.3~0.8 mmol/g	5g	29,200
	100~200 mesh, 1% DVB		25g	105,000
44901	H-Ile-2-Chlorotrityl Resin	0.3~0.8 mmol/g	5g	8,400
	100~200 mesh, 1% DVB		25g	33,900
44904	H-Ile-2-Chlorotrityl Resin	0.3~0.8 mmol/g	5g	8,400
	200~400 mesh, 1% DVB		25g	33,900
44902	Fmoc-Ile-2-Chlorotrityl Resin	0.3~0.8 mmol/g	25g	15,200
	100~200 mesh, 1% DVB		100g	52,500
44903	H-D-Allo-Ile-2-Chlorotrityl Resin	0.3~0.8 mmol/g	5g	29,200
	100~200 mesh, 1% DVB		25g	72,900
45001	H-Leu-2-Chlorotrityl Resin	0.3~0.8 mmol/g	5g	8,400
	100~200 mesh, 1% DVB		25g	33,900
45003	H-Leu-2-Chlorotrityl Resin	0.3~0.8 mmol/g	5g	8,400
	200~400 mesh, 1% DVB		25g	33,900
45006	H-Leucinol-2-Chlorotrityl Resin	0.3~0.8 mmol/g	5g	29,200
	200~400 mesh, 1% DVB		25g	105,000
45002	Fmoc-Leu-2-Chlorotrityl Resin	0.3~0.8 mmol/g	5g	8,400
	100~200 mesh, 1% DVB		25g	32,700
45005	H-D-Leu-2-Chlorotrityl Resin	0.3~0.8 mmol/g	5g	14,000
	100~200 mesh, 1% DVB		25g	49,000
45104	H-Lys-2-Chlorotrityl Resin	0.3~0.8 mmol/g	5g	9,600
	100~200 mesh, 1% DVB		25g	33,600
45101	H-Lys(Boc)-2-Chlorotrityl Resin	0.3~0.8 mmol/g	5g	9,600
	100~200 mesh, 1% DVB		25g	33,600
45107	H-Lys(Boc)-2-Chlorotrityl Resin	0.3~0.8 mmol/g	5g	9,600
	200~400 mesh, 1% DVB		25g	33,600
45109	H-Lys(Dde)-2-Chlorotrityl Resin	0.3~0.8 mmol/g	5g	58,300
	100~200 mesh, 1% DVB		25g	96,800
45103	Fmoc-Lys(Boc)-2-Chlorotrityl Resin	0.3~0.8 mmol/g	5g	9,300
	100~200 mesh, 1% DVB		25g	32,500
45108	Fmoc-Lys(Boc)-ol-2-Chlorotrityl Resin	0.3~0.8 mmol/g	5g	23,400
	100~200 mesh, 1% DVB		25g	70,000
45105	Fmoc-Lys(Dde)-2-Chlorotrityl Resin	0.3~0.8 mmol/g	5g	19,500
	100~200 mesh, 1% DVB		25g	58,300
45201	H-Met-2-Chlorotrityl Resin	0.3~0.8 mmol/g	5g	8,000
	100~200 mesh, 1% DVB		25g	32,000
45202	Fmoc-Met-2-Chlorotrityl Resin	0.3~0.8 mmol/g	5g	9,700
	100~200 mesh, 1% DVB		25g	31,500
45909	H-Orn(Boc)-2-Chlorotrityl Resin	0.3~0.8 mmol/g	5g	18,400
	100~200 mesh, 1% DVB		25g	64,200

Cat. #	品名	分子量	供給量	参考価格
45916	H-Orn(Boc)-2-Chlorotrityl Resin	0.3~0.8 mmol/g	5g	18,400
	200~400 mesh, 1% DVB		25g	64,200
45910	Fmoc-Orn(Boc)-2-Chlorotrityl Resin	0.3~0.8 mmol/g	5g	10,700
	100~200 mesh, 1% DVB		25g	37,200
45915	Fmoc-Orn(Boc)-2-Chlorotrityl Resin	0.3~0.8 mmol/g	5g	10,700
	200~400 mesh, 1% DVB		25g	37,200
45301	H-Phe-2-Chlorotrityl Resin	0.3~0.8 mmol/g	5g	8,400
	100~200 mesh, 1% DVB		25g	33,900
45307	H-Phe-2-Chlorotrityl Resin	0.3~0.8 mmol/g	5g	8,400
	200~400 mesh, 1% DVB		25g	33,900
45303	H-D-Phe-2-Chlorotrityl Resin	0.3~0.8 mmol/g	5g	8,700
	100~200 mesh, 1% DVB		25g	30,100
45304	Fmoc-N-Me-D-Phe-2-Chlorotrityl Resin	0.3~0.8 mmol/g	5g	11,700
	100~200 mesh, 1% DVB		25g	40,900
45911	Fmoc-D-Phg-2-Chlorotrityl Resin	0.3~0.8 mmol/g	5g	14,000
	100~200 mesh, 1% DVB		25g	46,700
45401	H-Pro-2-Chlorotrityl Resin	0.3~0.8 mmol/g	5g	8,000
	100~200 mesh, 1% DVB		25g	32,100
45406	Fmoc-Pro-2-Chlorotrityl Resin	0.3~0.8 mmol/g	5g	8,900
	100~200 mesh, 1% DVB		25g	30,900
45503	H-Ser(HPO3Bzl)-2-Chlorotrityl Resin	0.3~0.8 mmol/g	1g	6,900
	100~200 mesh, 1% DVB		5g	24,200
45004	Fmoc-Ser(tBu)-2-Chlorotrityl Resin	0.3~0.8 mmol/g	5g	9,400
	100~200 mesh, 1% DVB		25g	37,400
45505	Fmoc-Ser(Trt)-2-Chlorotrityl Resin	0.3~0.8 mmol/g	5g	12,300
	100~200 mesh, 1% DVB		25g	24,000
45501	H-Ser(tBu)-2-Chlorotrityl Resin	0.3~0.8 mmol/g	5g	9,400
	100~200 mesh, 1% DVB		25g	37,400
45502	H-Ser(Trt)-2-Chlorotrityl Resin	0.3~0.8 mmol/g	5g	12,900
	100~200 mesh, 1% DVB		25g	51,400
45506	H-D-Ser(tBu)-2-Chlorotrityl Resin	0.3~0.8 mmol/g	5g	29,200
	100~200 mesh, 1% DVB		25g	81,700
45919	H-Tic-2-Chlorotrityl Resin	0.3~0.8 mmol/g	5g	52,900
	100~200 mesh, 1% DVB		25g	88,100
45601	H-Thr(tBu)-2-Chlorotrityl Resin	0.3~0.8 mmol/g	5g	9,400
	100~200 mesh, 1% DVB		25g	37,400
45607	H-Thr(tBu)-2-Chlorotrityl Resin	0.3~0.8 mmol/g	5g	9,400
	200~400 mesh, 1% DVB		25g	37,400
45602	H-Thr(Trt)-2-Chlorotrityl Resin	0.3~0.8 mmol/g	5g	12,900
	100~200 mesh, 1% DVB		25g	51,400
45606	Fmoc-Thr(tBu)-2-Chlorotrityl resin	0.3~0.8 mmol/g	25g	17,500
	100~200 mesh, 1% DVB		100g	58,300
45608	Fmoc-Thr(tBu)-2-Chlorotrityl resin	0.3~0.8 mmol/g	5g	9,400
	200~400 mesh, 1% DVB		25g	37,400

Cat. #	品名	分子量	供給量	参考価格
45604	Fmoc-Thr(tBu)-ol-2-Chlorotrityl Resin	0.3~0.8 mmol/g	5g	29,200
	100~200 mesh, 1% DVB		25g	105,000
45701	H-Trp-2-Chlorotrityl Resin	0.3~0.8 mmol/g	5g	8,400
	100~200 mesh, 1% DVB		25g	38,500
45706	H-Trp-2-Chlorotrityl Resin	0.3~0.8 mmol/g	5g	8,400
	200~400 mesh, 1% DVB		25g	38,500
45702	H-Trp(Boc)-2-Chlorotrityl Resin	0.3~0.8 mmol/g	5g	10,000
	100~200 mesh, 1% DVB		25g	39,700
45703	Fmoc-Trp(Boc)-2-Chlorotrityl Resin	0.3~0.8 mmol/g	5g	10,000
	100~200 mesh, 1% DVB		25g	39,700
45804	Fmoc-Tyr(Bzl)-2-Chlorotrityl Resin	0.3~0.8 mmol/g	5g	12,900
	100~200 mesh, 1% DVB		25g	52,500
45803	Fmoc-Tyr(tBu)-2-Chlorotrityl Resin	0.3~0.8 mmol/g	5g	9,700
	100~200 mesh, 1% DVB		25g	29,200
45807	Fmoc-D-Tyr(Bzl)-2-Chlorotrityl Resin	0.3~0.8 mmol/g	5g	16,400
	100~200 mesh, 1% DVB		25g	56,000
45805	H-Tyr(Bzl)-2-Chlorotrityl Resin	0.3~0.8 mmol/g	5g	14,000
	100~200 mesh, 1% DVB		25g	54,900
45801	H-Tyr(tBu)-2-Chlorotrityl Resin	0.3~0.8 mmol/g	5g	9,500
	100~200 mesh, 1% DVB		25g	37,900
45806	H-D-Tyr(Bzl)-2-Chlorotrityl Resin	0.3~0.8 mmol/g	5g	17,500
	100~200 mesh, 1% DVB		25g	58,300
48107	H-D-Tyr(tBu)-2-Chlorotrityl Resin	0.3~0.8 mmol/g	5g	21,000
	100~200 mesh, 1% DVB		25g	73,500
45802	H-Tyr(Trt)-2-Chlorotrityl Resin	0.3~0.8 mmol/g	5g	12,900
	100~200 mesh, 1% DVB		25g	51,400
45912	Fmoc-Val-2-Chlorotrityl Resin	0.3~0.8 mmol/g	5g	8,200
	100~200 mesh, 1% DVB		25g	28,600
45901	H-Val-2-Chlorotrityl Resin	0.3~0.8 mmol/g	5g	8,400
	100~200 mesh, 1% DVB		25g	33,900
45917	H-Val-2-Chlorotrityl Resin	0.3~0.8 mmol/g	5g	8,400
	200~400 mesh, 1% DVB		25g	33,900