

TransBionovo provides a broad range of double-strand DNA molecular weight markers for conventional electrophoresis. All DNA markers are generated from restriction enzymes digested plasmids. All DNA markers are in ready-to-use format.

DNA Marker Selection Guide

DNA Marker	Agarose	DNA Marker	Agarose
<i>Trans2K</i> [®] DNA Marker	1.5%	1Kb DNA Ladder	1.0%
<i>Trans2K</i> [®] Plus DNA Marker	1.5%	1Kb Plus DNA Ladder	1.0%
<i>Trans2K</i> [®] Plus II DNA Marker	1.5%	100bp DNA Ladder	2.0%
<i>Trans5K</i> DNA Marker	1.5%	100bp Plus DNA Ladder	2.0%
<i>Trans15K</i> DNA Marker	0.8%	100bp Plus II DNA Ladder	2.0%

