

Item	Description	参考価格(税別)
APPLICATIONS KITS		
AppKit-01	NanoSpray Starter Kit: Packaged application kit including consumables needed for general nano-ESI analysis (platform independent)	¥579,000
AppKit-02	HT Profiling Application Kit-PicoChip Format	¥553,000
AppKit-03	Peptide Quantitation Application Kit	¥649,000
AppKit-04	Biomarker Discovery Application Kit..(TRC-25-Custom - (2) IT100-25H354)	¥598,000
AppKit-05	Top Down & PTM Application Kit..(TRC-25-Custom - (2) IT100-25H502)	¥800,000
AppKit-06	Metabolism Application Kit	¥856,000
AppKit-07	Metabolomics Application Kit	¥587,000
TUBING		
CT-360-25-SP	Pre-Cut Super Pack: Contains 25 pieces of 25 ID, Five pieces each of all five lengths	¥39,000
CT280-100-30-5	Precut Fused Silica Tubing: 280 OD, 100 ID, 30cm Length, 5 per pack	¥11,000
CT360-100-10-5	Precut Fused-Silica Tubing: 360 OD, 100 ID, 10cm Length(+/-1.0mm) 5 per pack	¥7,000
CT360-100-15-5	Precut Fused Silica Tubing: 360 OD, 100 ID, 15cm Length(+/-1.0mm) 5 per pack	¥8,000
CT360-100-25-5	Precut Fused Silica Tubing: 360 OD, 100 ID, 25cm Length(+/-1.0mm) 5 per pack	¥10,000
CT360-100-50-5	Precut Fused Silica Tubing, 360 OD, 100 ID, 50cm Length(+/-1.0mm) 5 per pack	¥14,000
CT360-100-6-5	Precut Fused Silica Tubing: 360 OD, 100 ID, 6cm Length(+/-1.0mm) 5 per pack	¥7,000
CT360-100-SP	Pre-Cut Super Pack: Contains 25 pieces of 100 ID, Five pieces each of all five lengths(+/-1.0mm)	¥37,000
CT360-20-10-5	Precut Fused Silica Tubing, 360 OD, 20 ID, 10cm Length(+/-1.0mm) 5 per pack	¥7,000
CT360-20-15-5	Precut Fused Silica Tubing, 360 OD, 20 ID, 15cm Length(+/-1.0mm) 5 per pack	¥8,000
CT360-20-25-5	Precut Fused Silica Tubing, 360 OD, 20 ID, 25cm Length(+/-1.0mm) 5 per pack	¥10,000
CT360-20-50-5	Precut Fused Silica Tubing, 360 OD, 20 ID, 50cm Length(+/-1.0mm) 5 per pack	¥14,000
CT360-20-6-5	Precut Fused Silica Tubing, 360 OD, 20 ID, 6cm Length(+/-1.0mm) 5 per pack	¥7,000
CT360-20-SP	Pre-Cut Super Pack: Contains 25 pieces of 20 ID, Five pieces each of all five lengths(+/-1.0mm)	¥37,000
CT360-50-10-5	Precut Fused Silica Tubing, 360 OD, 50 ID, 10cm Length(+/-1.0mm) 5 per pack	¥7,000
CT360-50-12-5	Pre-Cut Fused-Silica Tubing: 360 um OD, 50 um ID, 12 cm length(+/-1.0mm) 5 pieces per pack	¥8,000
CT360-50-15-5	Precut Fused Silica Tubing, 360 OD, 50 ID, 15cm Length(+/-1.0mm) 5 per pack	¥8,000
CT360-50-25-5	Precut Fused Silica Tubing, 360 OD, 50 ID, 25cm Length(+/-1.0mm) 5 per pack	¥10,000
CT360-50-45-5	Pre-cut fused-silica tubing: 360 um OD, 50 um ID, 45 cm length(+/-1.0mm) 5 pieces per pack	¥14,000
CT360-50-50-5	Precut Fused Silica Tubing, 360 OD, 50 ID, 50cm Length(+/-1.0mm) 5 per pack	¥14,000
CT360-50-6-5	Precut Fused Silica Tubing, 360 OD, 50 ID, 6cm Length(+/-1.0mm) 5 per pack	¥7,000
CT360-50-90-5	Pre-cut fused-silica tubing: 360 um OD, 50 um ID,90 cm length(+/-1.0mm), 5 pieces per pack	¥22,000
CT360-50-SP	Pre-Cut Super Pack: Contains 25 pieces of 50 ID, Five pieces each of all five lengths(+/-1.0mm)	¥37,000
CT360-75-10-5	Precut Fused Silica Tubing, 360 OD, 75 ID, 10cm Length(+/-1.0mm) 5 per pack	¥7,000
CT360-75-15-5	Precut Fused Silica Tubing, 360 OD, 75 ID, 15cm Length(+/-1.0mm) 5 per pack	¥8,000
CT360-75-25-5	Precut Fused Silica Tubing, 360 OD, 75 ID, 25cm Length(+/-1.0mm) 5 per pack	¥10,000
CT360-75-50-5	Precut Fused Silica Tubing, 360 OD, 75 ID, 50cm Length(+/-1.0mm) 5 per pack	¥14,000
CT360-75-6-5	Precut Fused Silica Tubing, 360 OD, 75 ID, 6cm Length(+/-1.0mm) 5 per pack	¥7,000
CT360-75-90-5	Pre-cut fused-silica tubing: 360 um OD, 75 um ID,90 cm length(+/-1.0mm) 5 pieces per pack	¥22,000
CT360-75-SP	Pre-Cut Super Pack: Contains 25 pieces of 75 ID, Five pieces each of all five lengths(+/-1.0mm)	¥37,000
RT240-100	Fused Silica Tubing: 240 OD, 100 ID, 2 meters	¥7,000
RT240-40	Fused Silica Tubing: 240 OD, 40 ID, 2 meters	¥7,000
RT285-100	Fused Silica Tubing: 285 OD, 100 ID, 2 meters	¥7,000
RT285-75	Fused Silica Tubing: 285 OD, 75 ID, 2 meters	¥7,000
RT360-10	Fused Silica Tubing: 360 OD, 10 ID, 2 meters	¥7,000
RT360-100	Fused Silica Tubing: 360 OD, 100 ID, 2 meters	¥7,000

価格・仕様は予告なく変更されている場合があります。

Item	Description	参考価格(税別)
RT360-100-1000	Fused Silica Tubing: 360 OD, 100 ID, 10 meters	¥30,000
RT360-150	Fused Silica Tubing: 360 OD, 150 ID, 2 meters	¥7,000
RT360-20	Fused Silica Tubing: 360 OD, 20 ID, 2 meters	¥7,000
RT360-20-1000	Fused Silica Tubing: 360 OD, 20 ID, 10 meters	¥30,000
RT360-200	Fused Silica Tubing: 360 OD, 200 ID, 2 meters	¥7,000
RT360-25	Fused Silica Tubing: 360 OD, 25 ID, 2 meters	¥7,000
RT360-50	Fused Silica Tubing: 360 OD, 50 ID, 2 meters	¥7,000
RT360-50-1000	Fused Silica Tubing: 360 OD, 50 ID, 10 meters	¥30,000
RT360-75	Fused Silica Tubing: 360 OD, 75 ID, 2 meters	¥7,000
RT360-75-1000	Fused Silica Tubing: 360 OD, 75 ID, 10 meters	¥30,000
RT360-KIT	Fused Silica Tubing: 2 meters of each , RT360-10, RT360-20, RT360-50, RT360-75, & RT360-100	¥30,000
TOOLS & ESSENTIALS		
FSFK-2	Fused Silica Fittings Kit with tubing assortment	¥70,000
MOUNT-FS1	Acrylic base for fused-silica inspection.	¥15,000
PS-STND-1	PicoSure Test Standard..Equimolar mixture of 7 synthetic peptides, 100 pmol	¥9,000
SCRIBE	Diamond wedge scribe	¥10,000
TIP-KIT	Nano-ESI Micro Tool Kit - Including - Fine forceps, medium forceps, diamond scribe, swaging needle and a 6" metal ruler	¥36,000
CONNECTORS & SLEEVES		
NO-C360ET	Tightening tool for 360µm OD tubing	¥2,000
NO-C360NFKPG	Custom PEEK nut/ferrule for FS 360µm OD tubing used w/Valco NO-C360US62PKG / NO-C360US66PKG	¥5,000
NO-C360US62PKG	Custom Valco stainless steel union for fused silica, 50um ID, w/PEEK nuts and ferrules	¥26,000
NO-C360US66PKG	Custom Valco stainless steel union for fused silica, 150um bore, Mw/PEEK nuts and ferrules	¥22,000
NO-NT1C18	NuTips Sample Prep Tips, C18	¥23,000
NO-ZU1TTI	ZDV (zero dead volume) Titanium Union; 1/16"	¥14,000
NO14-824	10 µL syringe with a 10 cm long, 170 µm OD flexible fused-silica needle. Ideal for filling offline PicoTips.	¥21,000
NO2-000-001	Calibrated Micropipets, 1-5 uL 250/box	¥10,000
NO22351656	Eppendorf GELoader Tips Box of 192 tips	¥27,000
NOC-1326	TRC Assembly - Nanobore Sample Trap Column Assembly with 1um Ti inlet frit (Upchurch C-1326) only (no columns)	¥148,000
NOC-NEU.5XFPK.35	Union, 100um bore, w/liners for 360 tubing	¥9,000
NOC-NL.35L-5	PEEK Liner (pack of 5)	¥3,000
NOC-NNFLFPK	Nanovolume Nuts and Ferrules for 10 Port Valve (package of 5)	¥10,000
NOC1100-01	Column holder part of trap cartridge assembly, PEEK, natural	¥12,000
NOC360RU.5FS2	SS reducing union, 1/32" to 360um, 50um thru bore	¥13,000
NOF-120BLK	1-Piece Fingertight fitting 1/16" ID 10-32 PEEK black	¥1,000
NOF-125x	Standard Head Nut 6-32 Coned, pkg of 10	¥12,000
NOF-142Nx	FingerTight ferrule, Tefzel (ETFE) 1/16" for 10/32 (10 pk)	¥8,000
NOF-152	Microferrule for 360 um OD tubing (Natural)	¥2,000
NOF-172x	Ferrule for Micro Tee/Cross, PEEK for 0.025" tubing, BLACK (pack of 10)	¥13,000
NOF-184x	MicroTight Sleeve, 330 um ID, ORANGE 10/pk	¥4,000
NOF-185x	MicroTight tubing sleeve, 395um ID for tubing size 340-380um, green, PEEK, 10 pk	¥4,000
NOF-192x	SealTight Ferrule, 10 pk	¥7,000
NOF-194x	Sealtight nut, 10-32 SS hexhead, includes F-192 PEEK ferrule for 1/16" OD tubing, 10pk	¥20,000
NOF-195x	SealTight nut / F-192 Peek Ferrule, 1/16", 10-32, short, peek black, 10 pk	¥11,000
NOF-231	Tubing sleeve 560µm ID x 1/16" OD PEEK natural	¥1,000
NOF-240x	FEP Teflon Sleeve, 0.011" ID x 1/16 OD, Blue 10/pack	¥4,000

Item	Description	参考価格(税別)
NOF-241x	SealTight (tm) FEP sleeve 1/16" OD, 300 um ID, orange, 10 per pack	¥4,000
NOF-242x	NanoTight tubing sleeve, FEP, green 395um ID x 1/16" OD for 340-380um tubing, pack of 10	¥4,000
NOF-330-01	Nut, fingertight, long 1/6", 10-32, PEEK (body only)	¥1,000
NOF-331Nx	NanoTight Nuts 10-32 w/142N Ferrule, 10 per pack	¥10,000
NOF-332-01	Nut, PEEK, natural	¥1,000
NOF-332Nx	SHORT HEX NANOTIGHT FITTINGW/F-142N FERRULE 10-32 10PK	¥12,000
NOF-345	WRENCH, FLUSHNUT, FOR 10-32 NUTS	¥15,000
NOF-350	Flushnut, SS,10-32 for 0.0625" tubing	¥2,000
NOF-381X	TUBING SLEEVE, PEEK, .007 IN (180µm) x 1/32 IN x 1.6 IN, YELLOW (10 PK)	¥4,000
NOF-385X	TUBING SLEEVE, PEEK, .015 IN (380µm) x 1/32 IN x 1.6 IN, NATURAL (10 PK)	¥4,000
NOM-120x	PEEK (tm) End-fittings for inline MicroFilter with .5 um frit, 10 per pack	¥21,000
NOM-124NF	Replacement capsule union, PEEK, no frit..For use with Sample Trap Cartridge Assy	¥4,000
NOM-126	Replacement Titanium Inlet Frit 2-pack (Upchurch M-126)	¥13,000
NOM-134	1 um NanoFilter Capsule, PEEK (tm) w/Titanium frit, 1/32", 2 per pack	¥13,000
NOM-436x	Ferrule Cluster for Micro Injection Valve, 10pk (Upchurch Valve model M-435)	¥10,000
NOM-436x	Replacemtn Ferrule Clusters for NanoPeak Valves (qty of 10)	¥11,000
NOM-520	In-line microfilter, 0.5 um, with (5) M-120 filters	¥17,000
NOM-538	Inline Nanofilter Assembly for 360 um tubing	¥45,000
NONT1C18	Nuttips C18 Tips 1-10 uL tips, Box of 96	¥27,000
NOP-200N	Flangeless ferrule 1/16" ID TEFZEL natural, used with NOP-201 on 450 sources.	¥400
NOP-200NX	Flangeless ferrule 1/16" ID TEFZEL natural, 10pk, used with NOP-201 on 450 sources.	¥2,000
NOP-201	1/16" Flangeless nut 1/4-28 Delrin, black; 1 piece	¥400
NOP-201X	1/16" Flangeless nut 1/4-28 Delrin, black, 10pk, used with NOP-200nx on the 450 source.	¥2,000
NOP-309	Plug 1/4-28 flat bottom, Delrin, black, 1 pc, used with the 450 source	¥1,000
NOP-309x	Plug 1/4-28 flat bottom, Delrin, black, 10 piece (PicoView 450)	¥3,000
NOP-629	Adapter, female slip Luer to female 10-32 TEFZEL, 1 pc	¥1,000
NOP-720	True ZDV union w/fittings and gauge plug	¥6,000
NOP-720-M1.0	Micro Tight Union, PEEK 6-32- Modified for 1.0mm Glass	¥9,000
NOP-720-M1.2	Micro Tight Union, PEEK 6-32- Modified for 1.2mm Glass	¥9,000
NOP-770	Upchurch P-770 union, PEEK, capillary to 1/6" OD	¥6,000
NOP-775	MicroTee, 150 um thru-holes, PEEK (tm)	¥19,000
NOP-777	MicroCross for MicroTight sleeves (150 um thru-holes)	¥22,000
NOP720M	Modified Micro Tight Union, PEEK 6-32	¥7,000
NOPV5510	Electronics housing assembly for DPV/PV-550 PicoView systems; sheath gas enabled	¥241,000
NOTNTC18	Trap'nTips: C18 sample prep tips for offline nanospray. Box of 96 tips.	¥53,000
NOTNTC4	Trap'nTips: C4 sample prep tips for offline nanospray. Box of 96 tip	¥62,000
NOTNTCAR	Trap'nTips: Carbon sample prep tips for offline nanospray. Box of 96 tips	¥62,000
NOTNTMC18	Trap'nTips: Carbon and C18 sample prep tips for offline nanospray. Box of 96 tips.	¥62,000
NOTNTTI	Trap'nTips: Titanium sample prep tips for offline nanospray. Box of 96 tips.	¥62,000
NOTNTZR	Trap'nTips: Zirconium sample prep tips for offline nanospray. Box of 96 tip	¥62,000
NOU-435-01	Union body, webbed, 10-32, 0.010" (.25mm) thru hole, SST	¥4,000
NOU-438	UNION, TRUE ZDV, WITH 2 U-400/401, 10-32, .067 IN (1.70mm) THRU HOLE, SST	¥4,000
NOV-472-11X	NanoFerrules PEEK, for NanoPeak Valves (qty of 10)	¥10,000
TIF-126	Replacement Titanium Inlet Frit 2-pack (Upchurch M-126)	¥9,000
PICOCLEAR CONNECTORS		

Item	Description	参考価格(税別)
PCU-360	PicoClear Union for 360 um OD tubing	¥25,000
PCU-360-I	For PicoClear Union (PCU-360): (2) Replacement cores for 360 um OD tubing with swaging tool	¥13,000
PCU-360-NF	PicoClear Union (PCU-360) replacement ferrules. Incl: (2) PEEK primary sealing ferrules, (2) PEEK secondary locking ferrules	¥11,000